

Product Information

PLA-PEG-PLA

Cat. No.: **S01YF-0723-KX11**

This product is for research use only and is not intended for diagnostic use.

Product Type

Phospholipids

Classification

PLA-PEG

Appearance

White solid

Molecular Weight

PLA MW: 1000 2000 5000 10000 15000 20000 40000; PEG MW: 1000 2000 3400 5000

Solubility

Chloroform; Dichloromethane

Purity

95%

Size

1g; 5g; 10g

Storage

-20°C

Storage Shelf Time

1 year