

---

# Product Information

## PLA-PEG-Iodoacetamide

Cat. No.: **S01YF-0723-KX38**

This product is for research use only and is not intended for diagnostic use.

### Product Type

Phospholipids

### Classification

PLA-PEG

### Appearance

White solid

### Molecular Weight

PLA MW: 1000 2000 5000 10000 15000 20000 40000; PEG MW: 1000 2000 3400 5000

### Solubility

Chloroform; Dichloromethane

### Purity

95%

### Size

500mg; 1g

### Storage

-20°C

### Storage Shelf Time

1 year