

Product Information

PLA-PEG-HA (Hyaluronate)

Cat. No.: **S01YF-0723-KX41**

This product is for research use only and is not intended for diagnostic use.

Product Type

Phospholipids

Classification

PLA-PEG

Appearance

White solid

Molecular Weight

PLA MW: 1000 2000 5000 10000 15000 20000 40000; PEG MW: 1000 2000 3400 5000

Solubility

Dichloromethane; Water

Purity

95%

Size

100mg; 500mg; 1g

Storage

-20°C

Storage Shelf Time

1 year