

Product Information

β-Cyclodextrin-PLA-PEG-PEI

Cat. No.: **S01YF-0723-KX49**

This product is for research use only and is not intended for diagnostic use.

Product Type

Phospholipids

Classification

PLA-PEG

Appearance

Solid

Molecular Weight

PEI MW: 600 1000 2000 5000; PEG MW: 1000 2000 3400 5000; PLA MW: 1000 2000 5000 10000

Purity

95%

Size

500mg

Storage

-20°C

Storage Shelf Time

1 year