

Product Information

PLA8k-TK-PEG2k-MA

Cat. No.: **S01YF-0723-KX57**

This product is for research use only and is not intended for diagnostic use.

Product Type

Phospholipids

Classification

PLA-PEG

Appearance

White solid

Molecular Weight

PEG MW: 350 550 750 1000 2000 3400 5000 10000 20000 30000 40000

Solubility

Chloroform; Methanol; DMSO

Purity

95%

Size

500mg; 1g

Storage

-20°C

Storage Shelf Time

1 year