

Product Information

PLA-PEG-ACUPA, MW: 2K, 2K

Cat. No.: **S01YF-0723-KX80**

This product is for research use only and is not intended for diagnostic use.

Product Type

Phospholipids

Classification

PLA-PEG

Appearance

White solid

Molecular Weight

PEG MW: 350 550 750 1000 2000 3400 5000 10000 20000 30000 40000

Solubility

Chloroform; Methanol; DMSO

Purity

95%

Size

1g

Storage

-20°C

Storage Shelf Time

1 year