

Product Information

PCL-PEG-Folate, MW: 3500, 2000

Cat. No.: **S01YF-0823-KX10**

This product is for research use only and is not intended for diagnostic use.

This product is a PCL-PEG phospholipid for research use only.

Product Type

Phospholipids

Classification

PCL-PEG

Appearance

Yellowish-brown solid

Molecular Weight

PCL3500, PEG2000

Solubility

Dichloromethane

Purity

Substitution=92%, PDI=1.10

Size

100mg; 500mg; 1g

Storage

-20°C

Storage Shelf Time

1 year