

Product Information

PCL-PEG-Multi-Arms (4 arm 8 arm)

Cat. No.: **S01YF-0823-KX113**

This product is for research use only and is not intended for diagnostic use.

This product is a PCL-PEG phospholipid for research use only.

Product Type

Phospholipids

Classification

PCL-PEG

Appearance

Solid

Molecular Weight

PCL: 1000 2000 5000 10000 15000 20000 40000; PEG: 1000 2000 3400 5000

Purity

>95%

Storage

-20°C in the dark

Storage Shelf Time

1 year