

Product Information

PCL-PEG-NH₂, MW:PCL2000/PEG2000

Cat. No.: **S01YF-0823-KX12**

This product is for research use only and is not intended for diagnostic use.

This product is a PCL-PEG phospholipid for research use only.

Product Type

Phospholipids

Classification

PCL-PEG

Appearance

White solid

Molecular Weight

PEG2000, PCL2000

Solubility

Dichloromethane

Purity

>95%

Size

100mg; 500mg; 1g

Storage

-20°C

Storage Shelf Time

1 year