

# Product Information

## PCL2000-Hyd-PEG2000-COOH

Cat. No.: **S01YF-0823-KX14**

This product is for research use only and is not intended for diagnostic use.

This product is a PCL-PEG phospholipid for research use only.

### Product Type

Phospholipids

### Classification

PCL-PEG

### Appearance

Yellow solid or paste

### Molecular Weight

PCL2000, PEG2000

### Solubility

Dichloromethane

### Purity

PDI<1.4

### Size

100mg; 500mg; 1g

### Storage

-20°C

### Storage Shelf Time

1 year