

Product Information

PCL-mPEG, MW: 10000, 2000

Cat. No.: **S01YF-0823-KX15**

This product is for research use only and is not intended for diagnostic use.

This product is a PCL-PEG phospholipid for research use only.

Product Type

Phospholipids

Classification

PCL-PEG

Appearance

White solid

Molecular Weight

PEG2000, PCL10000

Solubility

Chloroform; Dichloromethane

Purity

PDI<1.4

Size

1g; 5g; 10g

Storage

-20°C

Storage Shelf Time

1 year