

Product Information

mPEG-Hyd-PCL

Cat. No.: **S01YF-0823-KX19**

This product is for research use only and is not intended for diagnostic use.

This product is a PCL-PEG phospholipid for research use only.

Product Type

Phospholipids

Classification

PCL-PEG

Appearance

White solid

Molecular Weight

PEG: 1000 2000 5000 10000; PCL: 1000 2000 5000 10000 15000 20000 40000

Solubility

Dichloromethane

Purity

PDI<1.5

Size

500mg; 1g

Storage

-20°C

Storage Shelf Time

1 year