

# Product Information

## **PCL-PEG-Folate MW:3400,2000**

Cat. No.: **S01YF-0823-KX2**

This product is for research use only and is not intended for diagnostic use.

This product is a PCL-PEG phospholipid for research use only.

### **Product Type**

Phospholipids

### **Classification**

PCL-PEG

### **Appearance**

Yellowish-brown solid

### **Molecular Weight**

PCL3400, PEG2000

### **Solubility**

Chloroform; Dichloromethane

### **Purity**

Substitution>90%, PDI<1.4

### **Size**

100mg; 500mg; 1g

### **Storage**

-20°C

### **Storage Shelf Time**

1 year