

Product Information

PLGA-PEG-NH₂, MW:2000,3400, PLGA:50/50

Cat. No.: **S01YF-0823-KX243**

This product is for research use only and is not intended for diagnostic use.

This product is a PLGA-PEG phospholipid for research use only.

Product Type

Phospholipids

Classification

PLGA-PEG

Appearance

White solid

Molecular Weight

PEG3400, PLGA2000

Solubility

Chloroform; Dichloromethane

Purity

Substitution>90%, PDI<1.4

Size

100mg; 500mg; 1g

Storage

-20°C

Storage Shelf Time

1 year