

Product Information

mPEG-SS-PLGA-COOH MW:5000,5000,PLGA:50/50

Cat. No.: **S01YF-0823-KX244**

This product is for research use only and is not intended for diagnostic use.

This product is a PLGA-PEG phospholipid for research use only.

Product Type

Phospholipids

Classification

PLGA-PEG

Appearance

White solid

Molecular Weight

PEG5000, PLGA5000

Solubility

Chloroform; Dichloromethane

Purity

Substitution>90%, PDI<1.4

Size

200mg; 500mg

Storage

-20°C

Storage Shelf Time

1 year