

Product Information

PLGA-SH PLGA:50/50 MW:5000

Cat. No.: **S01YF-0823-KX245**

This product is for research use only and is not intended for diagnostic use.

This product is a PLGA-PEG phospholipid for research use only.

Product Type

Phospholipids

Classification

PLGA-PEG

Appearance

White solid

Molecular Weight

5100

Solubility

Chloroform; Dichloromethane

Purity

Substitution>90%, PDI<1.8

Size

500mg; 1g

Storage

-20°C

Storage Shelf Time

1 year