

Product Information

mPEG2000-S-S-PLGA1000-COOH PLGA:50/50

Cat. No.: **S01YF-0823-KX246**

This product is for research use only and is not intended for diagnostic use.

This product is a PLGA-PEG phospholipid for research use only.

Product Type

Phospholipids

Classification

PLGA-PEG

Appearance

Light yellow solid

Molecular Weight

PLGA1000, PEG2000

Solubility

Dichloromethane

Purity

Substitution=95%, PDI<1.4

Size

500mg; 1g

Storage

-20°C

Storage Shelf Time

1 year