

Product Information

PLGA-PEG-MAL,MW(PLGA45000,PEG5000) 50:50

Cat. No.: **S01YF-0823-KX247**

This product is for research use only and is not intended for diagnostic use.

This product is a PLGA-PEG phospholipid for research use only.

Product Type

Phospholipids

Classification

PLGA-PEG

Appearance

White solid

Molecular Weight

PLGA45000, PEG5000

Solubility

Dichloromethane

Purity

Substitution>90%, PDI<1.5

Size

100mg; 500mg; 1g

Storage

-20°C

Storage Shelf Time

1 year