

Product Information

PLGA-mPEG(MW,PLGA:45000 PEG:2000) 50:50

Cat. No.: **S01YF-0823-KX248**

This product is for research use only and is not intended for diagnostic use.

This product is a PLGA-PEG phospholipid for research use only.

Product Type

Phospholipids

Classification

PLGA-PEG

Appearance

White solid

Molecular Weight

PLGA45000, PEG2000

Solubility

Dichloromethane

Purity

PDI<1.4

Size

1 g; 5g

Storage

-20°C

Storage Shelf Time

1 year