

Product Information

PLGA-PEG-NHS, PLGA20000 PEG3400 50:50

Cat. No.: **S01YF-0823-KX249**

This product is for research use only and is not intended for diagnostic use.

This product is a PLGA-PEG phospholipid for research use only.

Product Type

Phospholipids

Classification

PLGA-PEG

Appearance

White solid

Molecular Weight

PLGA20000, PEG3400

Solubility

Dichloromethane

Purity

>90%

Size

100mg; 500mg; 1g

Storage

-20°C

Storage Shelf Time

1 year