

---

# Product Information

## **mPEG5000-b-PLGA10000**

Cat. No.: **S01YF-0823-KX256**

This product is for research use only and is not intended for diagnostic use.

This product is a PLGA-PEG phospholipid for research use only.

### **Product Type**

Phospholipids

### **Classification**

PLGA-PEG

### **Appearance**

White solid

### **Molecular Weight**

15000

### **Solubility**

Chloroform; Dichloromethane

### **Purity**

Substitution>90%, PDI<1.4

### **Size**

1g; 5g; 10g

### **Storage**

-20°C

### **Storage Shelf Time**

1 year