

# Product Information

## PCL-PEG-FITC

Cat. No.: **S01YF-0823-KX27**

This product is for research use only and is not intended for diagnostic use.

This product is a PCL-PEG phospholipid for research use only.

### Product Type

Phospholipids

### Classification

PCL-PEG

### Appearance

Green powder

### Molecular Weight

PEG1000 2000 3400 5000

### Solubility

Dichloromethane

### Purity

>95%

### Size

100mg; 500mg; 1g

### Storage

-20°C

### Storage Shelf Time

1 year