

Product Information

PLGA-PEG-Folate

Cat. No.: **S01YF-0823-KX275**

This product is for research use only and is not intended for diagnostic use.

This product is a PLGA-PEG phospholipid for research use only.

Product Type

Phospholipids

Classification

PLGA-PEG

Appearance

Yellow solid

Molecular Weight

PLGA: 1000, 2000, 5000, 10000, 15000; PEG: 2000, 3400, 5000

Solubility

Chloroform; Dichloromethane

Purity

Substitution>90%, PDI<1.4

Size

100mg; 500mg; 1g

Storage

-20°C

Storage Shelf Time

1 year