

Product Information

mPEG-b-PLGA

Cat. No.: **S01YF-0823-KX278**

This product is for research use only and is not intended for diagnostic use.

This product is a PLGA-PEG phospholipid for research use only.

Product Type

Phospholipids

Classification

PLGA-PEG

Appearance

White solid

Molecular Weight

22000; 4400; 8400; 5400; 23400

Solubility

Chloroform; Dichloromethane

Purity

Substitution>90%, PDI<1.4

Size

1g; 5g; 10g

Storage

-20°C

Storage Shelf Time

1 year