

# Product Information

## mPEG-SS-PLGA-NHS

Cat. No.: **S01YF-0823-KX294**

This product is for research use only and is not intended for diagnostic use.

This product is a PLGA-PEG phospholipid for research use only.

### Product Type

Phospholipids

### Classification

PLGA-PEG

### Appearance

White solid

### Molecular Weight

PLGA: 1000 2000 5000 10000 15000 20000 40000; PEG: 1000 2000 3400 5000

### Solubility

Chloroform; Dichloromethane

### Purity

Substitution>90%, PDI<1.4

### Size

100mg; 500mg; 1g

### Storage

-20°C

### Storage Shelf Time

1 year