

Product Information

mPEG-SS-PLGA-MAL

Cat. No.: **S01YF-0823-KX295**

This product is for research use only and is not intended for diagnostic use.

This product is a PLGA-PEG phospholipid for research use only.

Product Type

Phospholipids

Classification

PLGA-PEG

Appearance

White solid

Molecular Weight

PLGA: 1000 2000 5000 10000 15000 20000 40000; PEG: 1000 2000 3400 5000

Solubility

Chloroform; Dichloromethane

Purity

Substitution>90%, PDI<1.4

Size

100mg; 500mg; 1g

Storage

-20°C

Storage Shelf Time

1 year