

Product Information

PCL-Rhodamine B

Cat. No.: **S01YF-0823-KX31**

This product is for research use only and is not intended for diagnostic use.

This product is a PCL-PEG phospholipid for research use only.

Product Type

Phospholipids

Classification

PCL-PEG

Appearance

Red solid

Molecular Weight

PCL: 1000 2000 5000 10000 15000 20000 40000

Solubility

Dichloromethane

Purity

>95%

Size

100mg; 250mg; 500mg

Storage

-20°C

Storage Shelf Time

1 year