

Product Information

PLGA-PEG-Glycyrrhetic acid

Cat. No.: **S01YF-0823-KX317**

This product is for research use only and is not intended for diagnostic use.

This product is a PLGA-PEG phospholipid for research use only.

Product Type

Phospholipids

Classification

PLGA-PEG

Appearance

Solid

Molecular Weight

PLGA5000, PEG2000

Solubility

Dichloromethane

Purity

>95%

Size

100mg; 500mg; 1g

Storage

-20°C

Storage Shelf Time

1 year