

Product Information

PLGA-SS-PEG-AANL-KLVFFK(GGDGR)YIGSR MW:7600,2000 PLGA:75/25

Cat. No.: **S01YF-0823-KX342**

This product is for research use only and is not intended for diagnostic use.

This product is a PLGA-PEG phospholipid for research use only.

Product Type

Phospholipids

Classification

PLGA-PEG

Appearance

White solid

Molecular Weight

PEG: 350 550 750 1000 2000 3400 5000 10000 20000 30000 40000

Solubility

Chloroform; Methanol; DMSO

Purity

>95%

Size

1 g

Storage

-20°C

Storage Shelf Time

1 year