

Product Information

PLGA-PEG-Iodoacetate

Cat. No.: **S01YF-0823-KX362**

This product is for research use only and is not intended for diagnostic use.

This product is a PLGA-PEG phospholipid for research use only.

Product Type

Phospholipids

Classification

PLGA-PEG

Appearance

Solid

Molecular Weight

PLGA: 1000 2000 5000 10000 15000 20000 40000; PEG: 1000 2000 3400 5000

Purity

>95%

Storage

-20°C in the dark

Storage Shelf Time

1 year