

Product Information

PCL-PEG-OMe MW:2000,3400

Cat. No.: **S01YF-0823-KX4**

This product is for research use only and is not intended for diagnostic use.

This product is a PCL-PEG phospholipid for research use only.

Product Type

Phospholipids

Classification

PCL-PEG

Appearance

White solid

Molecular Weight

PCL3400, PEG2000

Solubility

Chloroform; Dichloromethane

Purity

Substitution>90%, PDI<1.4

Size

1g; 5g; 10g

Storage

-20°C

Storage Shelf Time

1 year