

Product Information

PCL-PEG-MAL, MW: PCL2000, PEG2000

Cat. No.: **S01YF-0823-KX5**

This product is for research use only and is not intended for diagnostic use.

This product is a PCL-PEG phospholipid for research use only.

Product Type

Phospholipids

Classification

PCL-PEG

Appearance

White solid

Molecular Weight

PCL2000, PEG2000

Solubility

Dichloromethane

Purity

Substitution>90%, PDI<1.3

Size

100mg; 500mg; 1g

Storage

-20°C

Storage Shelf Time

1 year