

Product Information

PCL-PEG-CO-MAL

Cat. No.: **S01YF-0823-KX71**

This product is for research use only and is not intended for diagnostic use.

This product is a PCL-PEG phospholipid for research use only.

Product Type

Phospholipids

Classification

PCL-PEG

Appearance

White solid

Molecular Weight

PCL: 1000 2000 5000 10000 15000 20000 40000; PEG: 1000 2000 3400 5000

Solubility

Chloroform; Dichloromethane

Purity

Substitution>90%, PDI<1.4

Size

200mg; 500mg; 1g

Storage

-20°C

Storage Shelf Time

1 year