

Product Information

PLL(20)-g[3.5]-PEG(2)

Cat. No.: **S01YF-0823-KX722**

This product is for research use only and is not intended for diagnostic use.

This product is a PLL-PEG phospholipid for research use only.

Product Type

Phospholipids

Classification

PLL-PEG

Appearance

Yellow solid

Molecular Weight

PLL20000, PEG2000

Solubility

Organic solvents

Purity

>90%

Size

100mg; 200mg; 500mg; 1g

Storage

-20°C

Storage Shelf Time

1 year