

Product Information

PLL2000-PEG5000, Poly(L-lysine)-b-PEG

Cat. No.: **S01YF-0823-KX726**

This product is for research use only and is not intended for diagnostic use.

This product is a PLL-PEG phospholipid for research use only.

Product Type

Phospholipids

Classification

PLL-PEG

Appearance

White solid

Molecular Weight

PLL2000, PEG5000

Solubility

Dichloromethane

Purity

>90%

Size

200mg; 500mg; 1g

Storage

-20°C

Storage Shelf Time

1 year