

Product Information

mPEG5000-PASP6000, Poly-L-aspartic acid-b-PEG

Cat. No.: **S01YF-0823-KX727**

This product is for research use only and is not intended for diagnostic use.

This product is a PLL-PEG phospholipid for research use only.

Product Type

Phospholipids

Classification

PLL-PEG

Appearance

White solid

Molecular Weight

PASP6000, PEG5000

Solubility

Dichloromethane

Purity

>90%

Size

200mg; 500mg; 1g

Storage

-20°C

Storage Shelf Time

1 year