

Product Information

PLL-PEG-NH2

Cat. No.: **S01YF-0823-KX729**

This product is for research use only and is not intended for diagnostic use.

This product is a PLL-PEG phospholipid for research use only.

Product Type

Phospholipids

Classification

PLL-PEG

Appearance

White solid

Molecular Weight

PLL: 1600 2500 4800 7400 12000 16000 24000 33000 48000, PEG: 1000 2000 3400 5000

Solubility

Chloroform; Dichloromethane

Purity

Substitution>90%, PDI<1.4

Size

200mg; 500mg

Storage

-20°C

Storage Shelf Time

1 year