

---

# Product Information

## mPEG-SS-PCL MW:5K,3K

Cat. No.: **S01YF-0823-KX76**

This product is for research use only and is not intended for diagnostic use.

This product is a PCL-PEG phospholipid for research use only.

### Product Type

Phospholipids

### Classification

PCL-PEG

### Appearance

White solid

### Molecular Weight

PEG: 350 550 750 1000 2000 3400 5000 10000 20000 30000 40000

### Solubility

Chloroform; Methanol; DMSO

### Purity

>95%

### Size

500 mg

### Storage

-20°C

### Storage Shelf Time

1 year