

Product Information

PEG-SS-PCL-PEI (Straight)

Cat. No.: **S01YF-0823-KX78**

This product is for research use only and is not intended for diagnostic use.

This product is a PCL-PEG phospholipid for research use only.

Product Type

Phospholipids

Classification

PCL-PEG

Appearance

Solid

Molecular Weight

PEG: 1000 2000 3400 5000; PEI: 600 1000 2000 10000; PCL: 1000 2000 5000 10000 15000 20000 40000

Purity

>95%

Size

200 mg

Storage

-20°C

Storage Shelf Time

1 year