

Product Information

PCL-mPEG, MW: 3500, 2000

Cat. No.: **S01YF-0823-KX8**

This product is for research use only and is not intended for diagnostic use.

This product is a PCL-PEG phospholipid for research use only.

Product Type

Phospholipids

Classification

PCL-PEG

Appearance

White solid

Molecular Weight

PCL3500, PEG2000

Solubility

Dichloromethane

Purity

PDI<1.2

Size

1 g; 2 g

Storage

-20°C

Storage Shelf Time

1 year