

Product Information

NHS-PEG2K-cRGDfK(RRRRRRRRRRC)

Cat. No.: **S01YF-0823-KX85**

This product is for research use only and is not intended for diagnostic use.

This product is a PCL-PEG phospholipid for research use only.

Product Type

Phospholipids

Classification

PCL-PEG

Appearance

White solid

Molecular Weight

PEG: 350 550 750 1000 2000 3400 5000 10000 20000 30000 40000

Solubility

Chloroform; Methanol; DMSO

Purity

>95%

Size

100 mg

Storage

-20°C

Storage Shelf Time

1 year