
Product Information

PCL-Dextran, MW:7500,4000

Cat. No.: **S01YF-0823-KX9**

This product is for research use only and is not intended for diagnostic use.

This product is a PCL-PEG phospholipid for research use only.

Product Type

Phospholipids

Classification

PCL-PEG

Appearance

White solid

Molecular Weight

PCL7500, Dextran4000

Solubility

Dichloromethane

Purity

Substitution>90%, PDI<1.2

Size

1 g; 2 g

Storage

-20°C

Storage Shelf Time

1 year