

---

# Product Information

## **PCL-PEOz-NH2 MW:2K,5K**

Cat. No.: **S01YF-0823-KX90**

This product is for research use only and is not intended for diagnostic use.

This product is a PCL-PEG phospholipid for research use only.

### **Product Type**

Phospholipids

### **Classification**

PCL-PEG

### **Appearance**

White solid

### **Solubility**

Chloroform; Methanol; DMSO

### **Purity**

>95%

### **Size**

100mg; 500mg; 1g

### **Storage**

-20°C