

---

# Product Information

## PCL-PEG-Acetylthio

Cat. No.: **S01YF-0823-KX97**

This product is for research use only and is not intended for diagnostic use.

This product is a PCL-PEG phospholipid for research use only.

### Product Type

Phospholipids

### Classification

PCL-PEG

### Appearance

Solid

### Molecular Weight

PCL: 1000 2000 5000 10000 15000 20000 40000; PEG: 1000 2000 3400 5000

### Purity

>95%

### Storage

-20°C in the dark

### Storage Shelf Time

1 year