

Product Information

Gold Nanoparticles, Au-mPEG5000

Cat. No.: **S01YF-0923-KX188**

This product is for research use only and is not intended for diagnostic use.

This product is a PEGylated Gold Nanoparticle for research use only.

Product Type

Phospholipids

Classification

PEGylated Gold Nanoparticle

Appearance

Red liquid

CAS Number

7440-57-5

Molecular Weight

Particle size: 15nm, 30nm, 50nm, 150nm

Solubility

Insoluble and dispersed in water.

Purity

Particle size deviation: ± 2 nm

Size

25 mL; 50 mL; 100 mL

Storage

2-8 °C

Storage Shelf Time

6 months