

Product Information

Gold Nanoparticles, Au-PEG550-COOH

Cat. No.: **S01YF-0923-KX189**

This product is for research use only and is not intended for diagnostic use.

This product is a PEGylated Gold Nanoparticle for research use only.

Product Type

Phospholipids

Classification

PEGylated Gold Nanoparticle

Appearance

Red liquid

CAS Number

7440-57-5

Molecular Weight

Particle size: 10nm, 20nm

Solubility

Insoluble and dispersed in water.

Purity

Particle size deviation: ± 2 nm

Size

25 mL; 50 mL

Storage

2-8 °C

Storage Shelf Time

6 months