

Product Information

Dextran modified Gold Nanoparticles (15nm)

Cat. No.: **S01YF-0923-KX287**

This product is for research use only and is not intended for diagnostic use.

This product is a Dextran Gold Nanoparticle for research use only.

Product Type

Phospholipids

Classification

Dextran Gold Nanoparticles

Appearance

Red liquid

CAS Number

7440-57-5

Molecular Weight

Dextran: 4k, 10k, 20k, 40k, 70k

Solubility

Insoluble and dispersed in water.

Purity

Particle size deviation: $\pm 2\text{nm}$

Size

25 mL; 50 mL

Storage

2-8 °C

Storage Shelf Time

6 months