

Product Information

AuNRs-PEG-NH₂, Aspect Ratio: 2.5

Cat. No.: **S01YF-0923-KX320**

This product is for research use only and is not intended for diagnostic use.

This product is a PEGylated Gold Nanorod for research use only.

Product Type

Phospholipids

Classification

PEGylated Gold Nanorods

Appearance

Blue liquid

CAS Number

7440-57-5

Molecular Weight

PEG: 550, 1k, 2k, 5k, 10k, 20k

Solubility

Insoluble and dispersed in water.

Size

25 mL

Storage

2-8 °C

Storage Shelf Time

6 months