

Product Information

Mesoporous silica coating Fe₃O₄ nanoparticles (50nm)

Cat. No.: **S01YF-0923-KX553**

This product is for research use only and is not intended for diagnostic use.

This product is a Mes silica Fe₃O₄ nano for research use only.

Product Type

Phospholipids

Classification

Mes silica Fe₃O₄ nano

Appearance

Black colloid

CAS Number

1317-61-9

Solubility

Insoluble in water

Purity

Particle size deviation: ±10nm

Size

5 mL; 10 mL

Storage

4-6 °C

Storage Shelf Time

6 months