

Product Information

MemDX™ Human CCR8 HEK293 Cell Line with cAMP-biosensor Reporter

Cat. No.: **S01YF-1222-KX873**

This product is for research use only and is not intended for diagnostic use.

Product Information

Target Protein

CCR8

Target Protein Species

Human

Accession Number

NM_005201.4

Protein Tag

Tag-free

Host Cell Type

HEK293

Target Classification

GPCR

Target Family

Chemokine Receptor

Target Research Area

Cancer Research; Infectious Research; Inflammation Research

Related Diseases

Molluscum Contagiosum; Kaposi Sarcoma

Product Properties

Morphology

Epithelial

Assay Types

Drug screening and biological assays

Assay Reporter

cAMP-biosensor

Resistance

Puromycin; Hygromycin

Stability

10 passages

Mycoplasma Testing

Negative

Biosafety Level

Level 1

Activity

Yes

Quantity

5x10⁶ cells

Form

Frozen cells

Freeze Medium

70% RPMI 1640 + 20% FBS + 10% DMSO

Culture Medium

DMEM + 10% FBS + 100µg/mL Hygromycin + 9µg/mL Puromycin

Selective Antibiotic(s)

Regular antibiotics active against mycoplasmas, bacteria and fungi.

Handling Notes

Frozen cells should be thawed immediately upon receipt and grown according to handling procedure to ensure cell viability and proper assay performance.

Note: Do not freeze the cells upon receipt as it may result in irreversible damage to the cell line.

Disclaimer: We cannot guarantee cell viability if the cells are not thawed immediately upon receipt and grown according to handling procedure.

Incubation

37°C with 5% CO₂

Applications

Drug screening and biological assays

Application Notes

Cells were plated in a 384-well plate and incubated overnight at 37°C and 5% CO₂ to allow the cells to attach and grow. Cells were then stimulated with a control for high-throughput drugs screening and functional assays.

Use Restrictions

These cells are distributed for research use only.

Shipping

Dry ice

Storage

Liquid nitrogen

Target

Full Name

C-C motif chemokine receptor 8

Introduction

This gene encodes a member of the beta chemokine receptor family, which is predicted to be a seven transmembrane protein similar to G protein-coupled receptors. Chemokines and their receptors are important for the migration of various cell types into the inflammatory sites. This receptor protein preferentially expresses in the thymus. I-309, thymus activation-regulated cytokine (TARC) and macrophage inflammatory protein-1 beta (MIP-1 beta) have been identified as ligands of this receptor. Studies of this receptor and its ligands suggested its role in regulation of monocyte chemotaxis and thymic cell apoptosis. More specifically, this receptor may contribute to the proper positioning of activated T cells within the antigenic challenge sites and specialized areas of lymphoid tissues. This gene is located at the chemokine receptor gene cluster region.

GPCR Signaling Pathway

The endogenous ligand is CC chemokines. Targeted protein activation can cause binding of Gi to Go protein which, in turn, cause an inhibition of adenylate cyclase and then decrease of cAMP concentration.

G coupling

Gi & Go

Endogenous Ligand

CC chemokines

Alternative Names

CY6; TER1; CCR-8; CKRL1; CDw198; CMKBR8; GPRCY6; CMKBRL2; CC-CKR-8; C-C chemokine receptor type 8; CC chemokine receptor 8; CC chemokine receptor CHEMR1; CC-chemokine receptor chemr1; chemokine (C-C motif) receptor 8; chemokine (C-C) receptor 8; chemokine (C-C) receptor-like 2; chemokine receptor-like 1; CCR8; C-C motif chemokine receptor 8

Gene ID

[1237](#)

UniProt ID

[P51685](#)